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EASYspin Plus Plant RNA Kit

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1. 桃果实、花、根、叶：Isolation, characterisation and phylogenetic analysis of resistance gene analogues in a wild species of peach (*Prunus kansuensis*).*Canadian Journal of Plant Science*, 2011, 91(6): 961-970
2. 樱桃花、叶、萼等各部位：Over-expression of the PaAP1 gene from sweet cherry (*Prunus avium* L.) causes early floweri.*Journal of Plant Physiology*, 2012, Available online 1 December 2012
3. 洋葱根、茎、蕾、叶、雌雄蕊等各部位：Cloning and Expression Analysis of A Putative B Class MADS-box Gene of AcPI in Onion. *Scientia Agricultura Sinica*, 2012, 45(23):4759-4769

4. 茺菁 : Isolation and Functional Characterisation of the Genes Encoding $\Delta 8$ -Sphingolipid Desaturase from *Brassica rapa*. *Journal of Genetics and Genomics* Volume 39, Issue 1, January 2012, Pages 47–59
5. 茺菁 1 : EXPRESSION, DIVERGENCE AND EVOLUTION OF THE CALEOSIN GENE FAMILY IN *BRASSICA RAPA*. *Arch. Biol. Sci.*, Belgrade, 65 (3), 863-876, 2013 DOI:10.2298/ABS1303863H
6. 番茄叶: Effect of Low Temperature Stress on the Expression of ProDH Gene and the Activities of the Proline Dehydrogenase in Leaves of Tomato Seedling. *Chinese Agricultural Science Bulletin* 2012,28(10):132-135
7. 梔子叶 : Isolation of High Quality Total RNA from *Gardenia jasminoides* Eills.*Chinese Agricultural Science Bulletin*.2012, 28(27):194-198
8. 油桐果实 : Cui Qinqin, Han Xiaojiao, Chen Yicun, Zhan Zhiyong, Lin Liyuan, Wang Yangdong. Isolation and Expression Characteristics of Biotin Carboxyl Carrier Protein Coding Gene (VfBCCP) from *Vernicia fordii*.*SCIENTIA SILVAE SINICAE*. 2012, 48(8): Available online August
9. 油桐果实 1 : Selection of Reliable Reference Genes for Gene Expression Studies Using Real-Time PCR in Tung Tree during Seed Development. *PLoS ONE*, 2012, 7(8): e43084
10. 紫菜 : Molecular cloning and expression analysis of ribosomal protein S7 gene from *Porphyra haitanensis*. *JOURNAL OF FISHERIES OF CHINA*, 2011, 35 (12) : 1814-1821
11. 石斛 : Molecular characterization of a mitogen-activated protein kinase gene DoMPK1 in *Dendrobium officinale*. *Acta Pharmaceutica Sinica*, 2012, 47 (12): 1703-1709
12. 石斛 1 : ESTs Analysis Reveals Putative Genes Involved in Symbiotic Seed Germination in *Dendrobium officinale*. *Symbiotic Germination Genes in D. officinale*. August 2013 | Volume 8 | Issue 8 | e72705
13. 大豆 : RNA-seq Analysis Reveals Ethylene-Mediated Reproductive Organ Development and Abscission in Soybean(*Glycine max* L. Merr.). *Plant Mol Biol Rep*, 2012, published online: 4 Dec, 2012
14. 大豆 1: Construction of ethylene regulatory network based on the phytohormones related gene transcriptome profiling and prediction of transcription factor activities in soybean. *Acta Physiologiae Plantarum*, 2012, published online: 12 Dec, 2012
15. 红花玉兰 : Expression Analysis of MAwAG in Different Organs and Developmental Stages of *Magnolia wufengensis*. *Chinese Bulletin of Botany*, 2013, 48 (2): 1–5
16. 毛桃 : Cloning and Phylogeny Analysis of PpAP2 Floral Homologous Genes in Peach. *Chinese Agricultural Science Bulletin*, 2013, 29(7): 99-104
17. 五倍子 : Cloning and characterisation of a phenylalanine ammonia-lyase gene from *Rhus chinensis*. *Plant Cell Rep*, 2013, published online: 15 March, 2013
18. :五倍子 1: Cloning, characterization and expression of chalcone synthase from medicinal plant *Rhus chinensis*.*J. Plant Biochem. Biotechnol.* DOI 10.1007/s13562-013-0231-9
19. 青杆 : cDNA Cloning and Bioinformatic Analysis of the sPPa1 Gene form *Picea wilsonii*. *Plant Science Journal*, 2012, 30(40): 394-401
20. 青杆 1: cDNA Cloning and Bioinformatic Analysis of PsbO Gene from *Picea wilsonii*.*Life Science Research*, 2012, 16(3): 201-206
21. 青杆 2: Cloning and Tissue Expression Analysis of PwPSAF in *Picea wilsonii*. *SCIENTIA SILVAE SINICAE*. Vol. 49, No. 10, Oct. 2013.

22. 洋葱: Molecular Cloning and Transcriptional Analysis of the Putative AGAMOUS Homolog AcAG in Onion (*Allium cepa*). *Plant Mol Biol Rep*, DOI 10.1007/s11105-013-0607-y
23. 木瓜: XsFAD2 gene encodes the enzyme responsible for the high linoleic acid content in oil accumulated in *Xanthoceras sorbifolia* seeds. *JOURNAL ARTICLE*. 2013-6-17.
24. 木瓜 1: Two novel diacylglycerol acyltransferase genes from *Xanthoceras 2 sorbifolia* are responsible for its seed oil content. *GENE-38688*; No. of pages: 9; 4C:
25. 柑橘: Efficient auto-excision of a selectable marker gene from transgenic citrus by combining the Cre/loxP system and ipt selection. *Plant Cell Rep*, DOI 10.1007/s00299-013-1470-x
26. 柑橘 1: Expression Analysis of Three Phloem-specific Promoters in Transgenic *Poncirus trifoliata*. *Acta Horticulturae Sinica*. 2014, 41(1): 1–8.
27. 柑橘 2: Activation of three pathogen-inducible promoters in transgenic citrus (*Citrus sinensis* Osbeck) after *Xanthomonas axonopodis* pv. citri infection and wounding. *Plant Cell Tiss Organ Cult*. DOI 10.1007/s11240-013-0423-y.
28. 茶梅花瓣: Comparison and Analysis of Methods of Extracting Total RNA from Petals of *Camellia sasanqua*. *Chinese Agricultural Science Bulletin*. 2013, 29(28):129-133.
29. 梔子: Isolation of High Quality Total RNA from *Gardenia jasminoides* Eills. *Chinese Agricultural Science Bulletin*. 2012, 28(27):194-198
30. 丹参: Genome-wide analysis and molecular dissection of the SPL gene family in *Salvia miltiorrhiza*. 2014 Jan;56(1):38-50. doi: 10.1111/jjpb.12111. Epub 2013 Nov 20.
31. 牡丹: Transcriptome Comparison Reveals Key Candidate Genes Responsible for the Unusual Reblooming Trait in Tree Peonies. *Genes Responsible for Reblooming in Tree Peonies*. November 2013 | Volume 8 | Issue 11 | e79996
32. 东南景天: Role of sulfur assimilation pathway in cadmium hyperaccumulation by *Sedum alfredii* Hance. *Ecotoxicology and Environmental Safety*. Volume 100, February 2014, Pages 159–165.
33. 山苍子: Identification of appropriate reference genes for normalizing transcript expression by quantitative real-time PCR in *Litsea cubeba*. *TECHNICAL NOTE*. *Mol Genet Genomics* (2013) 288:727–737, DOI 10.1007/s00438-013-0785-1
34. 木本植物: Heterologous gene silencing induced by tobacco rattle virus (TRV) is efficient for pursuing functional genomics studies in woody plants. *ORIGINAL PAPER*. *Plant Cell Tiss Organ Cult*, DOI 10.1007/s11240-013-0393-0
35. 棉花: Analysis of sea-island cotton and upland cotton in response to *Verticillium dahliae* infection by RNA sequencing. Sun et al. *BMC Genomics* 2013, 14:852 /1471-2164/14/852.
36. 桃子: Biochemical changes and defence responses during the development of peach gummosis caused by *Lasiodiplodia theobromae*. *Eur J Plant Pathol* (2014) 138:195–207, DOI 10.1007/s10658-013-0322-4.
37. 桃子 1: Carbohydrate metabolism changes in *Prunus persica* gummosis infected with *Lasiodiplodia theobromae*. *Phytopathology* "First Look" paper • <http://dx.doi.org/10.1094/PHYTO-01-13-0025-R> • posted 11/27/2013.
38. 海棠: The *Malus crabapple* transcription factor McMYB10 regulates anthocyanin biosynthesis during petal coloration. *Scientia Horticulturae* 166 (2014) 42–49.
39. 海藻: A rapid and sensitive method for field detection of *Prorocentrum donghaiense* using reverse transcription-coupled loop-mediated isothermal amplification. *Harmful Algae* 29 (2013) 31–39.